



December 14th 2022

OpenViBE: an open source BCI
software suite

PART 2 - Designing BCIs with OpenViBE

Arthur Desbois, Marie-Constance Corsi

ARAMIS team, Paris Brain Institute

PART 2 - Designing BCI systems with OpenViBE


2.1 - The Motor Imagery Paradigm

- **Underlying idea:**


Taking advantage of a **neurophysiological phenomenon** to establish a communication between the brain and the computer

Illustration with **Motor imagery-based BCI**


Motor execution VS Rest



Motor imagery VS Rest




Power decrease








Desynchronization effect
(Pfurtscheller et al, 1999)



- Behavioral properties
 - Movement / preparation for movement : Event-related desynchronization (ERD) (Pfurtscheller, G, Lopes da Silva, FH, 1999)
 - With relaxation/post-movement period : ERS
- Why using it in BCI ?
 - Mu/Beta activity modulation by motor-imagery, a way to communicate
 - Use of power spectra
 - To establish this communication :
 - Spatial selection
 - Frequency selection




Illustrations from BCI2000 website

Motor Imagery - In Practice...












PART 2 - Designing BCI systems with OpenViBE

2.2 - Technical Demo






Step 1 - Warmup! Simple I/O & Display





1. Search for “generic stream reader” in boxes list
2. Drag & drop to designer window
3. Set filename with browser
`C:\openvibe-3.4.0-64bit\share\openvibe\scenarios\signals\bci-motor-imagery.ov`

Step 1 - Warmup! Simple I/O & Display




4. Search for “signal display” in boxes list
5. Drag & drop to designer window
6. Drag & drop connection between boxes, using the “signal” stream type
7. Press Play!




Step 2 - DSP, Visualization

- Features to extract (recap)
 - Power Spectra
 - Sensorimotor Area
 - Mu/Alpha (8-12Hz) &/or Beta (14-29Hz)




(Gonzalez-Astudillo et al, 2020)

(Maeder et al., 2012)



- **OpenViBE** can manage multiple signal streams in parallel
 - In our example, 1 channel = 1 EEG electrode
 - We will load a pre-recorded signal file, that used 11 labeled electrodes and consider only a sub-set of those electrodes.

Step 2 - DSP, Visualization




Note:
Right click + Rename
can be useful when
designing a scenario...

Display windows will use
« renamed » label

1. Import a Generic Stream Reader Box, set the filename to:
`<install folder>\share\openvibe\scenarios\signals\bci-motor-imagery.ov`
2. Import a Channel Selector Box, link the boxes together
3. Set the selection to:
`C3;C4;FC3;FC4;C5;C1;C2;C6;CP3;CP4`
4. Display the output

- We now need to compute power spectra for the alpha & beta bands, for the selected electrodes

Step 2 - DSP, Visualization



Configure Temporal Filter settings

Filter Method:

Filter Type:

Filter Order:

Low Cut-off Frequency (Hz):

High Cut-off Frequency (Hz):

Band Pass Ripple (dB):

Override settings with configuration file

Configure Time based epoching settings

Epoch duration (in sec):




Epoch intervals (in sec):

Override settings with configuration file


1. “Temporal Filter”, set to [8;12]Hz
2. “Time based epoching”, of 1s, every 0.2s (overlapping windows of signal for power computation)
3. DSP #1
formula $x*x$
4. Average
(averaging x^2 across a window of 1s)
5. DSP #2
formula $\log_{10}(1+x)$

 $\Rightarrow \log_{10}(1+\text{avg}(x^2))$
6. Add Displays and rename them to “alpha band power” and “alpha band signal”


Step 2 - DSP, Visualization




1. Use the “widget reordering tool” for ease of use!
2. Set time scales : 10 (seconds) for alpha band signal display, 50 for alpha band power (since we have 5 times more data to display)
3. Set vertical scaling to “global”




Step 2 - DSP, Visualization




Step 2 - DSP, Visualization




- Repeat for Beta Band [12;29hz]
- Don't hesitate to copy/paste boxes...
... but make sure you check/update their parameters



Step 2 - DSP, Visualization



Alternative boxes for spectral estimation & visualization...




- **“Autoregressive Coefficients”**

Estimates an AR model for the signal

Output #1 = the AR coefficients



Output #2 = the PSD of the signal, estimated as the FFT of the AR model.



- **“Instant Bars”**

Allows for quickly visualizing spectra / vectors / etc...



⚠ a bit bulky to setup



- **Goal reminder:**

We want to train a classification algorithm...

- ...to distinguish between mental states (MI / REST) ...
- ... using EEG signals from the sensorimotor areas ...
- ... more precisely instantaneous spectral power...
- ... in bands alpha and/or beta



(Gonzalez-Astudillo et al, 2020)

(Maeder et al., 2012)

- **Goal reminder:**

We want to train a classification algorithm...

- ...to distinguish between mental states (MI / REST) ...
- ... using EEG signals from the sensorimotor areas ...
- ... more precisely instantaneous spectral power...
- ... in bands alpha and/or beta

- **In OpenViBE:**

Step 3 - Training the classifier

- Goal reminder:

We want to train a classification algorithm...

- ...to distinguish between mental states (MI / REST) ...
- ... using EEG signals from the sensorimotor areas ...
- ... more precisely instantaneous spectral power...
- ... in bands alpha and/or beta

- In OpenViBE:

- Feature Aggregator
- Classifier Trainer

- **Goal reminder:**

We want to train a classification algorithm...

- ...to distinguish between mental states (MI / REST) ...
- ... using EEG signals from the sensorimotor areas ...
- ... more precisely instantaneous spectral power...
- ... in bands alpha and/or beta

- **In OpenViBE:**

- Stimulation Based Epoching
- Feature Aggregator
- Classifier Trainer

- **Goal reminder:**

We want to train a classification algorithm...

- ...to distinguish between mental states (MI / REST) ...
- ... using EEG signals from the sensorimotor areas ...
- ... more precisely instantaneous spectral power...
- ... in bands alpha and/or beta

- **In OpenViBE:**

- Stimulation Based Epoching
- Channel Selector
- Feature Aggregator
- Classifier Trainer



- **Goal reminder:**

We want to train a classification algorithm...

- ...to distinguish between mental states (MI / REST) ...
- ... using EEG signals from the sensorimotor areas ...
- ... more precisely instantaneous spectral power...
- ... in bands alpha and/or beta

- **In OpenViBE:**

- Stimulation Based Epoching
- Channel Selector
- Time Based Epoching
- AutoRegressive Coefficients
- Feature Aggregator
- Classifier Trainer



- **Goal reminder:**

We want to train a classification algorithm...

- ...to distinguish between mental states (MI / REST) ...
- ... using EEG signals from the sensorimotor areas ...
- ... more precisely instantaneous spectral power...
- ... in bands alpha and/or beta




- **In OpenViBE:**

- Stimulation Based Epoching
- Channel Selector
- Time Based Epoching
- AutoRegressive Coefficients
- Frequency Band Selector
- Feature Aggregator
- Classifier Trainer


... + a few other boxes

Step 3 - Training the classifier



Load a signal file...

Step 3 - Training the classifier



Load a signal file...

- ... distinguish between mental states (MI / REST)...

Configure Stimulation based epoching settings


Epoch duration (in sec) ↑ ↓

Epoch offset (in sec) ↑ ↓

Stimulation to epoch from

☐ Override settings with configuration file

Step 3 - Training the classifier




Load a signal file...

• ... distinguish between mental states (MI / REST)...

• ... using EEG signals from the sensorimotor areas...

Step 3 - Training the classifier





Load a signal file...


... distinguish between mental states (MI / REST)...

... using EEG signals from the sensorimotor areas...

... more precisely instantaneous spectral power...




Step 3 - Training the classifier




Load a signal file...

- ... distinguish between mental states (MI / REST)...
- ... using EEG signals from the sensorimotor areas...
- ... more precisely instantaneous spectral power...
- ... in bands alpha and/or beta



(+ average across the freq. band)

Step 3 - Training the classifier



Load a signal file...

• ... distinguish between mental states (MI / REST)...

• ... using EEG signals from the sensorimotor areas...

• ... more precisely instantaneous spectral power...

• ... in bands alpha and/or beta

• ... to train the classification algo

Step 3 - Training the classifier

Configure Classifier trainer settings

Train trigger: OVTk_StimulationId_Train

Filename to save configuration to: \${Player_ScenarioDirectory}/motor-imagery-b

Multiclass strategy to apply: Native

Class 1 label: OVTk_GDF_Left

Class 2 label: OVTk_GDF_Right

Algorithm to use: Linear Discriminant Analysis (LDA)

Use shrinkage: false

Shrinkage coefficient (-1 == auto): -1.000000

Shrinkage: Force diagonal cov (DDA): false

Number of partitions for k-fold cross-validation test: 7


Balance classes: false

Override settings with configuration file

Buttons: Load..., Save..., Default, Revert, Apply, Cancel

- “Train Trigger”:
Set in the acquisition LUA script (see Annex) to a specific stimulation code, triggered at the end of the experiment
- Set the **path/filename** for the training “weights” (results)
- Check that the **classifier classes** are correctly labeled
- Select a **classification algo.** (LDA, SVM...) and set its parameters


Step 3 - Training the classifier



Other (useful!) boxes

- **Player Controller / Timeout:**
Control the flow of the experiment
- **Common Avg Reference**
Make sure the avg. of all channels' signals is zero
- **Epoch Average**
Average "Feature Vectors" used for classification across a full "Epoch" (here, the duration of a MI trial)
- **Identity**
Simply propagates signals & stims. Useful for readability



Step 3 - Training the classifier




```
Cross-validation test accuracy is 90.000000% (sigma = 20.000000%)
Cls vs cls      1      2
Target  1: 100.0   0.0 %, 10 examples
Target  2:  20.0  80.0 %, 10 examples
Training set accuracy is 100.000000% (optimistic)
Cls vs cls      1      2
Target  1: 100.0   0.0 %, 10 examples
Target  2:   0.0 100.0 %, 10 examples
```

Filename to save configuration to

output_weights.xml



Motor Imagery Protocol - Phase 2




- **Goals:**

From EEG data acquired in **real-time**, while providing **MI tasks** to the subject...

- We want to **classify their mental states** (also in real-time!)
- ... using the classifier previously trained.

➔ We will reuse most of the building blocks set up so far

Step 3 - Online classification




Using a pre-recorded signal file

In a real BCI experiment, we should of course use the **Acquisition Client** (and record the incoming signals for future replays and studies).

Today, we're running an online demonstration without EEG hardware, so we shall use a pre-recorded signal file...

The subject would receive instructions synchronized with stimulations set up in the **LUA stimulator Box**. Let's set this up exactly as for a real-life experiment.

Step 3 - Online classification




Signal Processing / Feature extraction

We want the classifier to receive the same type of data used for training.

Here, we shall use the same features previously used for training the classifier.

Step 3 - Online classification




Use the “**Classifier Processor**” box, loading the “**weights**” file generated in the training phase


Use the same **LUA script** & **Graz Visualization** boxes as for the **Acquisition phase** to generate stimulations & “tasks”, and display instructions for the subject.


The “**streamed matrix**” stream from **Classifier processor** to **Graz Visualization** corresponds to the “**classification accuracy**”

ACQUISITION SERVER/CLIENT PROTOCOL MANAGEMENT / STIMULATIONS




- **Scenario:**
in the workshop github:
`sc1-monitor-acq.xml`
- **OpenViBE Acquisition server:**
in the install folder:
`openvibe-acquisition-server.cmd`






1. In the “Driver” list, select “Generic Oscillator”
In this list, you’ll find all the drivers for OpenViBE’s supported EEG hardware
2. Click “Connect”
3. Click “Play”




EEG data acquisition


Stimulation generation
& display




- **LUA Stimulator Box**
(LUA = scripting language)
- **Experiment example:**
Different stimulation/event codes
at different times.
Useful for signal segmentation
("epoching")
- **Stimulation label examples:**
 - Experiment Start/Stop
 - Trial Start/Stop
 - LEFT / RIGHT
 - Button pressed
 - etc



- “Graz Visualization”
(specific for “Graz Protocol”,
based on Box “Display Cue
Image”)



- Displays specified image upon
receiving specified stimulation
code
- Transmits the **stimulation code &
time** to the Acquisition Server
- The stimulations are received
by the ACQ Client, synchronized
with the signal




- **“Player Controller”:**
Orchestrates the experiment course, by applying an action upon receiving a stimulation

GOING FURTHER...


SHARED PARAMETERS, METABOXES

Shared parameters / “Scenario configuration”




When using more complex scenarios, with parameters shared by multiple boxes...

- ... using **shared parameters** can ease manipulations and avoid making mistakes!



When using more complex scenarios, with parameters shared by multiple boxes...



- ... using **shared parameters** can ease manipulations and avoid making mistakes!
- **Set up shared parameters** in the “Scenario configuration” tab, next to “Boxes” in the right-part of the Designer GUI.
- Add/edit your parameters by clicking on “**Configure settings**”

Shared parameters / “Scenario configuration”

The screenshot displays a software interface for scenario configuration. On the left, a flowchart shows a 'Channel Selector CP3 In|Out|Set' box connected to four 'Time based epoching' boxes. A dialog box titled 'Configure Time based epoching settings' is open, showing 'Epoch duration (in sec)' set to '\$var{TimeEpochLength}' and 'Epoch intervals (in sec)' set to '\$var{TimeEpochShift}'. Below these fields are buttons for 'Load...', 'Save...', 'Default', 'Revert', 'Apply', and 'Cancel'. On the right, the 'Scenario Configuration' panel shows three settings: 'TimeEpochLength' (0.5), 'TimeEpochShift' (0.4), and 'ArModelOrder' (19). Each setting has up/down arrows and a copy icon.

Configure Time based epoching settings

Epoch duration (in sec)

Epoch intervals (in sec)


Override settings with configuration file

Scenario Configuration


TimeEpochLength


TimeEpochShift

ArModelOrder



- “Metaboxes” can be useful to replace recurring/repeating processing chains.
- Use them in conjunction with scenario configurations to simplify your scenarios







- Saving your metaboxes:


```
<install folder>\share\openvibe\metaboxes\
```

... they will appear in the “Metabox” folder the Box Browser next time you start the Designer.

- Full tutorial:

<http://openvibe.inria.fr/designer-tutorial-5-metaboxes/>





Configure Feature Extraction - Power Spectrum settings

ChannelSelected FC1

EpochDuration 1 ↑ ↓

EpochInterval 1 ↑ ↓

FilterOrder 10 ↑ ↓

PsdSize 0 ↑ ↓

FreqBand 15:22

EpochAvg 18 ↑ ↓

☒ Override settings with configuration file

Load... Save... Default Revert Apply Cancel

The metabox's parameters are the ones set in the “scenario configuration” of the metabox scenario

Boxes Scenario Configuration Scenario I/O

ChannelSelected FC1

EpochDuration 1 ↑ ↓

EpochInterval 1 ↑ ↓

FilterOrder 10 ↑ ↓

PsdSize 0 ↑ ↓

FreqBand 15:22

EpochAvg 18 ↑ ↓

Channel Selector In/Out/Set

Time based epoching



AutoRegressive Coefficients

Frequency Band Selector

Epoch average In/Out/Set

Feature aggregator In/Out/Set

TOPOGRAPHY VISUALIZATION



- Add “Electrode localization file reader” box, using file:

```
<install folder>\share\openvibe\  
electrode_sets\electrode_set_standard_cartesian.txt
```

- Add “2D Topographic map”
- Connect boxes, play the scenario

